

# AcuityWare Image Print 2

Installation Documentation

Feb 17, 2011



Copyright © 2011, ACUITYWARE, LLC. All rights reserved. No part of this documentation may be reproduced in any form or by any means or used to make any derivative work (such as translation, transformation, or adaptation) without written permission from ACUITYWARE, LLC.

ACUITYWARE, LLC reserves the right to revise this documentation and to make changes in content from time to time without obligation on the part of ACUITYWARE, LLC to provide notification of such revision or change.

ACUITYWARE, LLC provides this documentation without warranty of any kind, either implied or expressed, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. ACUITYWARE may make improvements or changes in the product(s) and/or the program(s) described in this documentation at any time.

All technical data and computer software is commercial in nature and developed solely at private expense. As the User, or Installer/Administrator of this software, you agree not to remove or deface any portion of any legend provided on any licensed program or documentation contained in, or delivered to you in conjunction with, this User Guide.


AcuityWare Imaging Package and the AcuityWare Image logo are trademarks of ACUITYWARE, LLC

Other brand and product names may be registered trademarks or trademarks of their respective holders.


Document version: 2.00


## Installation

To begin installation of the AcuityWare Image Print product, please double click on the executable that was delivered. You will see the following startup screen as depicted below. Go ahead and click “Next >” to proceed.




Please review the license agreement and click on “I accept the agreement” to continue installing this product.







The following screen prompts the installer for their Name and Company information. Please enter it into the fields provided and click "Next >" to continue.




The operator can then select the destination folder for the installation of the product. Make the appropriate changes and then proceed by clicking “Next >” to continue.




Review the final screen to make sure all items are correct. Once you have reviewed all items, click on the “Install” button to start the actual installation.



This screen will appear as the installation progresses. You can cancel the installation at any time during this process by clicking on the “Cancel” button.



Upon completion, the following screen will appear and prompt the installer to click the “Finish” button.



## Troubleshooting

### *FAQ (Frequently Asked Questions)*

- Q. Does AcuityWare image Print work with Windows 7?
- A. It has been repackaged and reworked to run in the Windows 7 environment.
- Q. Is it a 64 bit software package.
- A. No. It is designed to be run in the Windows 32 bit operating system. You can however try to run it on a 64 bit version and you might have to test it thoroughly prior to implementing it across the enterprise.
- Q. Why does the installation tell me that it already exists and I have to uninstall it before installing the package?
- A. The Imaging software was designed to be installed fresh. If you have an issue, uninstall it first before attempting to reinstall it. This way it will register the components and the software correctly within the operating system.

## Technical Support

If you are an existing AcuityWare Image Software customer and are having difficulties installing or using the software, you can find help using the following methods.

- Email: [support@acuityware.com](mailto:support@acuityware.com)
- Phone: 609-607-8114
- Fax: 609-994-3127

